

The family of fully integrated systems For Hyperthermic Intraperitoneal Chemotherapy

The ThermoChem™ System by ThermaSolutions was the first fully integrated system specifically designed and manufactured for Hyperthermic Intraperitoneal Chemotherapy (HIPEC). By exploiting heat and harnessing its power, this advanced medical technology offers a new choice in adjunctive surgical therapies. Used intraoperatively, the ThermoChem™ HT-2000 or HT-1000 Systems raise the temperature in the peritoneum by continuously circulating sterile solution with chemotherapy throughout the peritoneal cavity.



HIPEC Disposable Kit

The HIPEC Disposable kit is a complete HIPEC kit for a single procedure. The disposable is compatible with both the ThermoChem™ HT-1000 & HT-2000 devices. The single-use, sterile kit contains is pre-assembled and includes: Inflow & outflow cannula, temperature probes, heat transfer exchanger, flow reservoir with filtration, pressure sensor, waste bag and tubing & connectors to complete the circuit to the patient.

Order numbers:

TSP 1016:	1 inflow + 3 outflow
TSP 1019.0:	3 inflow + 2 outflow
TSP 2011:	2 inflow + 3 outflow

ThermaSolutions Europe BV

Takkebijsters 41
4817 BL Breda
THE NETHERLANDS

www.thermasolutions.com

Phone: +31 76 579 1144
Fax: +31 76 579 1149