

The family of fully integrated systems for Intraoperative Hyperthermia

The ThermoChem™ System by ThermaSolutions was the first fully integrated system specifically designed and manufactured for intraoperative hyperthermia. By exploiting heat and harnessing its power, this advanced medical technology offers a new choice in adjunctive surgical therapies. Used intraoperatively, the ThermoChem™ HT-2000 or HT-1000 Systems raise the temperature in the peritoneum by continuously circulating sterile solution throughout the peritoneal cavity.



IPH Disposable Kit

The IPH Disposable kit is a complete IPH kit for a single procedure. The disposable is compatible with both the ThermoChem™ HT-1000 & HT-2000 devices. The single-use, sterile kit contains pre-assembled and includes: Inflow & outflow cannulae, temperature probes, heat transfer exchanger, flow reservoir with filtration, pressure sensor, and tubing & Y-connectors to complete the circuit to the patient.

Order number # 111-2000

ThermaSolutions
1889 Buerkle Road
White Bear Lake, MN 55110

www.thermasolutions.com
Phone: +1 651 209 3900
Fax: +1 651 209 3903