


Isolated Limb Infusion Kit



Isolated Limb Infusion (ILI) Procedure Kit order number # ILI-8000

The ILI Procedure Kit is a complete single-use, sterile kit that contains pre-assembled and connectable components. The kit includes: inflow & outflow tubing, heat transfer exchanger, waste collection bags, three-way manifold to complete the circuit to the patient for regionally concentrated treatment.

ThermaSolutions
1889 Buerkle Road
White Bear Lake, MN 55110

www.thermasolutions.com
Phone: +1 651 209 3900
Fax: +1 651 209 3903