

ThermoCasing Tubing Wrap



ThermoCasing Tubing Wrap

Order number # ANC-3610
Contents: 3' single-use tubing wrap

The ThermoCasing product is a 1/2" elastomeric insulative tubing wrap for inflow and outflow tubing. Used between the IPH device and the sterile field, the ThermoCasing wrap surrounds both sets of tubing in operating room environments where ambient temperature may effect fluid temperature.

ThermaSolutions LLC
1889 Buerkle Road
White Bear Lake, MN 55110

www.thermasolutions.com
Phone: +1 651 209 3900
Fax: +1 651 209 3903