

Isolated Limb Infusion Kit

The ILI Procedure Kit is a complete single-use, sterile kit that contains pre-assembled and connectable components



The kit includes: inflow & outflow tubing, heat transfer exchanger, waste collection bags, three-way manifold to complete the circuit to the patient for regionally concentrated treatment

Isolated Limb Infusion (ILI) Procedure Kit
order number ILI-8000

ThermaSolutions
1889 Buerkle Road
White Bear Lake, MN 55110

www.thermasolutions.com
Phone: +1 651 209 3900
Fax: +1 651 209 3903