

# Isolated Limb Perfusion Kit



## Isolated Limb Perfusion (ILP) Procedure Kit order number # ILP-9001/ILP-9002

The ILP Procedure Kit is a complete single-use, sterile kit that contains is pre-assembled and connectable components. The kit includes: oxygenator, inflow & outflow tubing, heat transfer exchanger, flow reservoir with filtration, tubing & Y-connectors to complete the circuit to the patient.

**ThermaSolutions**  
1889 Buerkle Road  
White Bear Lake, MN 55110

**[www.thermasolutions.com](http://www.thermasolutions.com)**  
Phone: +1 651 209 3900  
Fax: +1 651 209 3903