


## Fully integrated system for Hyperthermic Intraperitoneal Chemotherapy (HIPEC)

The ThermoChem™ System by ThermoSolutions was the first fully integrated system specifically designed and manufactured for intraperitoneal hyperthermia. By exploiting heat and harnessing its power, this advanced medical technology offers a new choice in adjunctive surgical therapies. Used intraoperatively, the ThermoChem™ HT-2000 device and associated disposables, raises the temperature in the peritoneum by continuously circulating fluids throughout the peritoneal cavity as part of the HIPEC procedure.



### Lasso Tubing

The Lasso Tubing is a 3/8" outflow tubing that can be used with the ThermoSolutions IPH Disposable Kits (IPH-2000) and is compatible with both the ThermoChem™ HT-2000 and HT-1000 Systems. The Lasso Tubing acts as the outflow cannula, allowing flow of fluid from the patient back to the IPH circuit.

**ThermoSolutions**  
1889 Buerkle Road  
White Bear Lake, MN 55110

**Order number # FLW-0010**

[www.thermasolutions.com](http://www.thermasolutions.com)

Phone: +1 651 209 3900

Fax: +1 651 209 3903