


Fully integrated system for Hyperthermic Intraperitoneal Chemotherapy (HIPEC)

The ThermoChem™ System by ThermoSolutions was the first fully integrated system specifically designed and manufactured for intraperitoneal hyperthermia. By exploiting heat and harnessing its power, this advanced medical technology offers a new choice in adjunctive surgical therapies. Used intraoperatively, the ThermoChem™ HT-2000 device and associated disposables, raises the temperature in the peritoneum by continuously circulating fluids throughout the peritoneal cavity as part of the HIPEC procedure.



Long Lasso Tubing

The Long Lasso Tubing is a 3/8" outflow cannula that can be used with ThermoSolutions IPH Disposable Kits (IPH-2000) and is compatible with both the ThermoChem™ HT-2000 and HT-1000 Systems. The Long Lasso Tubing acts as the outflow cannula, allowing flow of fluid from the patient back to the IPH circuit. The tubing allows for easy placement in a minimally invasive case.

Order number # FLW-0111

ThermoSolutions
1889 Buerkle Road
White Bear Lake, MN 55110

www.thermasolutions.com
Phone: +1 651 209 3900
Fax: +1 651 209 3903