

## The family of fully integrated systems For Hyperthermic Intraperitoneal Chemotherapy

The ThermoChem™ Unit by ThermoSolutions was the first fully integrated system specifically designed and manufactured for Hyperthermic Intraperitoneal Chemotherapy (HIPEC). By exploiting heat and harnessing its power, this advanced medical technology offers a new choice in adjunctive surgical therapies. Used intraoperatively, the ThermoChem™ HT-2000 or HT-1000 Unit raise the temperature in the peritoneum by continuously circulating sterile solution with chemotherapy throughout the peritoneal cavity.



### HIPEC Hyperthermic Perfusion Kit

The HIPEC Hyperthermic Perfusion kit is a complete HIPEC kit for a single procedure. The disposable is compatible with both the ThermoChem™ HT-1000 & HT-2000 devices. The single-use, sterile kit contains pre-assembled and includes: Inflow & outflow cannula, temperature probes, heat transfer exchanger, flow reservoir with filtration, pressure sensor, waste bag and tubing & connectors to complete the circuit to the patient.

**Order numbers:** TSP 1016: 1 inflow + 3 outflow  
TSP 1019.0 2 inflow + 3 outflow drains outflow  
TSP 2011: 2 inflow + 3 outflow

**ThermoSolutions Europe BV**

Takkebijsters 41  
4817 BL Breda  
THE NETHERLANDS

**[www.thermasolutions.com](http://www.thermasolutions.com)**

Phone: +31 76 579 1144  
Fax: +31 76 579 1149